

# PREempt®

Concentrated One-Step  
Surface Cleaner and  
Disinfectant formulated

with **ACCELERATED®**  
HYDROGEN PEROXIDE  
Patented



Concentrate

## WHY CHOOSE PREempt™?

- **CLEANER** In-use solution is a one-step disinfectant cleaner. It has powerful surfactants, which means it disinfects while it cleans surfaces and equipment, resulting in lower cost and faster results.

---

- **FASTER** Realistic contact times on a broad spectrum of organisms, resulting in user compliance and confidence that disinfection has occurred.

---

- **RESPONSIBLE** Users not exposed to VOC (volatile organic compounds). In-use solution proven to be essentially non-toxic and non-irritating to eyes and skin.

---

- **SUSTAINABLE** Hydrogen Peroxide breaks down into water and oxygen, leaving no active residual and reducing environmental impact.

## DIRECTIONS FOR USE

Use Dilutions Made Easy

RATIO	PER 1 LITRE OF WATER	PER 4 LITRES OF WATER
1:40	25mL	100mL
1:128	8mL	32mL

**Disinfection of non-porous hard surfaces** coming in contact with intact skin such as the exterior of fume hoods, tabletops and other environmental surfaces.

1. Dilute 1:40.
2. Apply to surface with cloth or disposable wipe.
3. Ensure surface remains wet for 5 minutes.

*For use in laboratories, clean rooms and other critical environments that require cleaning and surface disinfection. This includes work stations, fume hoods, laboratory counter tops, equipment and other hard non-porous environmental surfaces.*



## AVAILABLE PRODUCTS

PREempt™ Concentrated One-Step Surface Cleaner and Disinfectant

**DIN:** 02436701

**SHELF LIFE:** 24 Months, 30 Days Diluted

**CONCENTRATION OF ACTIVE:** 7% Hydrogen Peroxide (Concentrated)

1. **Format Size:** 4 L Bottle  
**Product Code:** 11305  
**Packaging Specifications:** 4x4 L
2. **Format Size:** 200 L Drum  
**Product Code:** 11311  
**Packaging Specifications:** 1x200 L

## Disinfection of Non-Critical, Pre-cleaned Instruments and Medical Devices:

1. Dilute to 1:40.
2. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and dried.
3. Thoroughly clean and rinse lumens of hollow instruments.
4. Using either a tray or an ultrasonic unit, immerse instrument in use dilution for 5 minutes.
5. Remove and rinse instruments.
6. Wipe dry prior to use.
7. Discard solution after each use.

## EFFECTIVE AGAINST A BROAD-SPECTRUM OF BACTERIA, VIRUSES AND FUNGI IN 5 MINUTES!

### Virucidal: 5 Min. (1:40 dilution)

- Poliovirus Type 1, Sabin strain type 1 which allows for a Broad-Spectrum Virucide claim against most enveloped and non-enveloped viruses
- HIV-1 Human Immunodeficiency virus (HIV)
- Human Coronavirus 229E
- Feline Calicivirus, F9 Strain, as a surrogate for Norwalk and Norwalk-like viruses
- Effective against Poliovirus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus

### Bactericidal: 5 Min. (1:40 dilution)

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Vancomycin-resistant *Enterococcus faecalis* (VRE)

### Fungicidal: 5 Min. (1:40 dilution)

- *Trichophyton mentagrophytes*

### Tuberculocidal: 5 Min. (1:40 dilution)

- *Mycobacterium terrae*

### Broad-Spectrum Sanitizing: 30 Sec. (1:128 dilution)

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Klebsiella pneumoniae*
- *Escherichia coli* O157:H7
- *Campylobacter jejuni*



*PREempt™ Concentrate One-Step Surface Cleaner and Disinfectant is designed for use in many establishments including compound pharmacy labs, laboratories, clean rooms and critical environments where broad-spectrum germicidal performance in short contact time is required. PREempt™ Concentrate is specifically formulated for cleaning and disinfection of hard non-porous surfaces.*

*\*Please refer to reference sheet for a complete list of pathogenic organisms, additional features, and claims.*

*PREempt™ brand products are manufactured by Virox Technologies Inc.*

*Accelerated Hydrogen Peroxide® and PREempt™ are registered trademarks of Virox Technologies Inc.*



For more information on this or our other products, visit:

**virox.com**

2770 Coventry Road, Oakville, ON CANADA L6H 6R1  
1-800-387-7578 | [www.virox.com](http://www.virox.com)



*Engineering Revolutionary Disinfectants for the War Against Microbes*