

PREempt[®] plus

Wipes One-Step³ Disinfectant Cleaner

For Use in Health Care, Dental Offices Clinics, Pharmacies, and Laboratories

This product is a one-step³ germicidal disinfectant cleaner and deodorant designed for general cleaning and disinfecting of hard, non-porous surfaces.

It cleans quickly by removing dirt, grime, and other common soils and stains found in clinics, dental offices, pharmacies and laboratories.

This product is designed for use on hard, non-porous surfaces: aluminum, chrome, glazed ceramic, glazed porcelain, laminated surfaces, painted surfaces, plastic surfaces, stainless steel, and vinyl, associated with high touch surfaces, such as: carts, chairs, counters, diagnostic equipment, exam tables, or any hard, non-porous washable surface where disinfection is required.

CLAIMS

Bactericidal in 1 Minute:

Acinetobacter baumannii (Multi-Drug Resistant (MDR))
Bordetella bronchiseptica
Bordetella pertussis
Burkholderia cepacia
Enterobacter aerogenes
Enterobacter cloacae (New Delhi metallo-beta-lactamase 1 (NDM -1))
Enterococcus faecalis (Vancomycin Resistant (VRE))
Enterococcus faecium
Escherichia coli
Escherichia coli (extended spectrum beta-lactamase (ESBL))
Escherichia coli (New Delhi metallo-beta-lactamase 1 (NDM-1))
Escherichia coli O157:H7
Klebsiella pneumoniae (Carbapenem Resistant (CRE))
Klebsiella pneumoniae (extended spectrum beta-lactamase (ESBL))
Klebsiella pneumoniae (New Delhi metallo-beta-lactamase 1 (NDM-1))
Klebsiella pneumoniae
Legionella pneumoniae
Listeria monocytogenes
Pasteurella multocida
Pseudomonas aeruginosa
Pseudomonas aeruginosa (metallo-beta-lactamase (MBL) positive)
Salmonella enterica (formerly known as *Salmonella choleraesuis*)
Serratia marcescens
Shigella dysenteriae
Staphylococcus aureus
Staphylococcus aureus (Community Acquired Methicillin Resistant Genotype USA 400 (CA-MRSA 400))
Staphylococcus aureus (Methicillin Resistant (MRSA))
Staphylococcus aureus (Vancomycin Intermediate Resistant (VISA))
Staphylococcus aureus (Vancomycin Resistant (VRSA))
Staphylococcus epidermidis (Methicillin Resistant (MRSE))
Staphylococcus pseudintermedius
Stenotrophomonas (Pseudomonas) maltophilia
Streptococcus pneumoniae (Penicillin Resistant (PRSP))
Streptococcus pyogenes

Reference Sheet

KEEP OUT OF REACH
OF CHILDREN

Manufactured by:
Virox Technologies Inc.
2770 Coventry Road
Oakville, ON L6H 6R1
1-800-387-7578 | www.virox.com



USER SHOULD READ S.D.S.
BEFORE USING PRODUCT
SDS is available online at www.virox.com

Tuberculocidal in 1 Minute:

Mycobacterium bovis (BCG)

*Virucidal in 1 minute unless otherwise noted:

Adenovirus Type 8
Avian Influenza A Virus (H3N2)
Avian Influenza A Virus (H7N9)
Cytomegalovirus
Feline Calicivirus
Hepatitis B Virus (HBV) – 30 seconds
Hepatitis C Virus (HCV) – 30 seconds
Herpes Simplex Virus Type 1
Herpes Simplex Virus Type 2
HIV Type 1 Strain HTLV III_B – 30 seconds
Human Coronavirus
Influenza A Virus (Strain A/8/68)
Influenza B Virus
Norovirus
Pandemic Influenza A Virus H1N1
Parainfluenza Virus Type 3
Poliovirus Type 1 (Strain Chat)
Respiratory Syncytial Virus (RSV)
Rhinovirus Type-37
Rotavirus (Strain WA)
SARS-Related Coronavirus 2 (Human Coronavirus SARS-CoV-2
(the virus that causes COVID-19*))
Vaccinia Virus

Fungicidal in 1 minute:

Candida albicans
Microsporum canis
Trichophyton interdigitale (formerly known as *Trichophyton mentagrophytes*) – 30 seconds

Non-Food Contact Sanitizer in 10 Seconds against:

Acinetobacter baumannii
Campylobacter jejuni
Enterobacter aerogenes
(Vancomycin Resistant) *Enterococcus faecalis* (VRE)
Escherichia coli O157:H7
Listeria monocytogenes
Pseudomonas aeruginosa
Salmonella enterica (formerly known as *Salmonella choleraesuis*)
Staphylococcus aureus
(Methicillin Resistant) *Staphylococcus aureus* (MRSA)

Mold/Mildew (Mildewcidal Activity) at a 4 minute contact time:

Aspergillus brasiliensis (formerly known as *Aspergillus niger*)

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. A potable water rinse is required for food contact surfaces.

Not for personal cleansing. This is not a baby wipe. Do not flush in toilet.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the blood stream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

For Use as a One-Step Cleaner/Disinfectant:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous surfaces.
3. Allow surface to remain visibly wet for 1 minute.
4. Allow to air dry.

For Use as a Non-Food Contact Sanitizer:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain visibly wet for 10 seconds.
4. Allow to air dry.

To Kill Fungi:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow surface to remain visibly wet for 1 minute.
4. Allow to air dry.

For Use as a Tuberculocide:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow the surface to remain visibly wet for 1 minute.
4. Allow to air dry.

To Kill Mold and Mildew:

1. Pre-clean surfaces.
2. Pull towelette from canister and wipe product onto hard, non-porous surfaces until completely wet.
3. Allow the surface to remain visibly wet for 4 minutes.
4. Allow to air dry.

To Disinfect Non-Critical[‡], Pre-cleaned Instruments:

1. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and rough dried.
2. Wipe all hard, non-porous surfaces of instruments with this product until thoroughly wet.
3. Allow surfaces to remain visibly wet for 1 minute.
4. Allow to air dry.

To Kill *Trichophyton interdigitale*:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
3. Allow surface to remain visibly wet for 30 seconds.
4. Allow to air dry.

For Use to Clean and Disinfect Life Science Laboratory Surfaces, Instruments, and Glassware:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous surfaces.
3. Allow surface to remain visibly wet for 1 minute.
4. For glassware/instrument: Wipe surface thoroughly and let air dry before reuse.
For surfaces: Wipe surfaces dry or rinse.

To Disinfect Surfaces in Pharmaceutical and Cosmetic Processing Facilities:

1. Pre-clean visibly soiled areas.
2. Pull towelette from canister and wipe hard, non-porous surfaces.
3. Allow surface to remain visibly wet for 1 minute.
4. Rinse with sterile water.
5. Allow to air dry.

***KILLS HIV-1 AND HCV AND HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings (Hospitals, Nursing Homes, Dental Clinics) and other settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV OR HCV ON SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS:

Personal Protection: Disposable latex or vinyl gloves, gowns, face masks, and eye coverings must be worn during all cleaning of body fluids, blood, and decontamination procedures.

Cleaning Procedures: Blood/body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Allow surface to remain visibly wet for 30 seconds to kill HIV-1, HBV and HCV.

Disposal of Infectious Material: Cleaning materials used that may contain blood/body fluids must be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

Meets surface disinfection recommendations under the OSHA Bloodborne Pathogen Standard for *HBV, *HCV & *HIV-1

DISCLAIMER: You must read and follow the instructions on the label for our product and use the product for the particular applications specified on the product label. To the full extent permitted under applicable law, Virox Technologies Inc. shall not be liable for any damages whatsoever, arising out of or otherwise in connection with any failure to follow those instructions or any use of its products which is not specified on the product's label. The instructions on the product label are not a warranty, express or implied, and Virox Technologies, Inc. expressly disclaims any implied warranties, including any warranty of fitness for a particular purpose. Under applicable law, some of the above waivers or disclaimers of warranty, or exclusions of liability for certain types of damages may not apply to you. In such cases, our liability will be limited to the fullest extent permitted by applicable law.

[‡] Accelerated Hydrogen Peroxide® technology and designs are trademarks of Diversey, Inc. "Accelerated" does not imply increased product efficacy.

[±] Hard, non-porous surfaces

[‡] Non-Critical medical devices are items that come in contact only with intact skin.

³ When used according to disinfection use directions.