

PREempt[®]

One-Step Surface Cleaner
and Disinfectant Wipes
formulated with



Wipes

WHY CHOOSE PREempt™?

- **CLEANER** PREempt™ Wipes are a one-step disinfectant cleaner. It has powerful surfactants, which means it disinfects while it cleans surfaces and equipment, resulting in lower cost and faster results.
- **FASTER** Realistic contact times on a broad spectrum of organisms resulting in user compliance and confidence that disinfection has occurred.
- **RESPONSIBLE** Users are not exposed to VOC (volatile organic compounds). Proven to be essentially non-toxic and non-irritating to eyes and skin.
- **SUSTAINABLE** Hydrogen Peroxide breaks down into water and oxygen, leaving no active residual and reducing environmental impact.

DIRECTIONS FOR USE

Disinfection of non-porous hard surfaces coming in contact with intact skin such as the exterior of fume hoods, tabletops and other environmental surfaces.

1. Heavily soiled surfaces require cleaning prior to disinfection.
2. Apply to surface with disposable wipe.
3. Ensure surface remains wet for 3 minutes.

Disinfection of Non-Critical, Pre-cleaned Instruments and Medical Devices:

1. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and dried.
2. Wipe the instrument with PREempt™ Wipes and allow surface to remain wet for 3 minutes.
3. Wipe with a clean, damp cloth or paper towel and allow to air dry.

For use in laboratories, clean rooms and other critical environments that require cleaning and surface disinfection. This includes work stations, fume hoods, laboratory counter tops, equipment and other high touch non-porous environmental surfaces.

AVAILABLE PRODUCTS

PREempt™ One-Step Surface Cleaner and Disinfectant Wipes

DIN: 02436671

SHELF LIFE: 24 Months

CONCENTRATION OF ACTIVE:
0.5% Hydrogen Peroxide

1. **Format Size:** 6"x7" – 160ct Wipes
Product Code: 11221
Packaging Specifications: 12x1 Case
2. **Format Size:** 7.5"x10" – 10ct Wipes (Soft Pack)
Product Code: 11223
Packaging Specifications: 1 Case x 6 inner Cartons x 20 Packs

EFFECTIVE AGAINST A BROAD-SPECTRUM OF BACTERIA, VIRUSES AND FUNGI IN 3 MINUTES!

Virucidal: 3 Min.

- Poliovirus Type 1, Sabin strain type 1 which allows for a Broad-Spectrum Virucide claim against most enveloped and non-enveloped viruses
- HIV-1 Human Immunodeficiency virus (HIV)
- Human Coronavirus 229E
- Feline Calicivirus, F9 Strain, as a surrogate for Norwalk and Norwalk-like viruses
- Effective against Poliovirus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus

Bactericidal: 3 Min.

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Vancomycin-resistant *Enterococcus faecalis* (VRE)

Fungicidal: 3 Min.

- *Trichophyton mentagrophytes*

Tuberculocidal: 3 Min.

- *Mycobacterium terrae*

Broad-Spectrum Sanitizing: 30 Sec.

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Klebsiella pneumoniae*
- *Escherichia coli* O157:H7
- *Campylobacter jejuni*



PREempt™ One-Step Surface Cleaner and Disinfectant Wipes are designed for use in many establishments including compound pharmacy labs, laboratories, clean rooms and critical environments where the requirement for broad-spectrum germicidal performance in short contact time is required.

*Please refer to reference sheet for a complete list of pathogenic organisms, additional features, and claims.

PREempt™ brand products are manufactured by Virox Technologies Inc.

Accelerated Hydrogen Peroxide® and PREempt™ are registered trademarks of Virox Technologies Inc.



For more information on this or our other products, visit:

> virox.com

2770 Coventry Road, Oakville, ON CANADA L6H 6R1
1-800-387-7578 | www.virox.com



Engineering Revolutionary Disinfectants for the War Against Microbes